

ASEPTIX **DISINFECTION** SYSTEMS

GERM FREE

in a **FLASH**

C Diff

Viruses

MRSA

VRE

CLASS1**INC.**

Building Better Healthcare™

A COMBINED APPROACH FOR PERSISTENT,


ASEPTIX¹

Fixed UV Disinfection

FEATURES

- Automatically activates after each room use (when door closes)
- Infrared motion detectors prevent unwanted exposure
- Reactivates every 4 hours when room is not in use
- The most effective cleaning possible in **5 minutes**
- Tested to show **>99.9999% reduction** of *C. difficile*, MRSA, VRE and more

World's first fully automatic UVC room disinfection system
Ideal for patient bathrooms and equipment storage rooms



Faster than
any other
mobile UV
disinfection

ASEPTIX²

Mobile UV Disinfection

FEATURES

- Dual unit operation guarantees the most effective cleaning possible in **5 to 15 minutes**
- Minimizes “shadow area”
- Tested to show **>99.9999% reduction** of C. difficile, MRSA, VRE and more
- Integrated infrared motion detectors prevents unwanted exposure
- Data capture records disinfection process
- Wireless communication and tracking with any smart device

World's fastest mobile UVC room disinfection system
Ideal for terminal cleaning of patient rooms and ORs

Light weight dual system allows you to easily move ASEPTIX² between rooms



ASEPTIX² Disinfection Performance Data



	Distance	5'	6'	7'	8'	9'	10'
	mW/cm ²	1314	999	778	620	504	417
Klebsiella pneumoniae	Seconds	18	24	30	42	48	60
VRE	Seconds	30	36	42	54	66	78
MRSA	Minutes	1.2	1.6	2.1	2.6	3.2	3.8
C Diff Spores	Minutes	4.6	6.0	7.7	9.7	11.9	14.4

Log 6, 5 minutes

CLASS1 INC.

Building Better Healthcare™

Class 1 Inc.
 1 (800) 242 9723
info@class1inc.com
www.class1inc.com

Proud to be a founding member



CHAIR

Coalition for Healthcare
 Acquired Infection Reduction

chaircanada.org